

# Ocean Nutrient Analyzer

San José State University  
College of Engineering  
April 6, 2022

SJSU



Contributors: A. Silva, C. Hayes, T. Hunter, T. Ly





# Introduction and Presentation Overview

We are presenting updates on the following:

- Thermal Analysis
- Gantt Chart
- Port control
- Pump control
- Sequence control
- Reagent housing design
- Final enclosure design
- Data collection

# Thermal Analysis Update: We're going to be ok!

## Given:

AC-DC Power Supply Rating: 108 [W]

Stepper Motor Wattage Rating: 48 [W]

$P_{total} = 3 \cdot 48 + 108 = 348$  [W]

Power\_wall\_socket = 1800 [W]

**Find:** Temperature rise in enclosure

## Assumptions:

- 10% of wattage dissipates as heat

## Wattage Dissipated as Heat

$1800W - 348W = 1416$  [W]

$1416 \cdot .10 = 141.6$  Watts dissipated as heat

## Surface Area of Outside of Enclosure

$107.8125 \text{ in}^2 \cdot 2 = 215.625$  [in<sup>2</sup>]


$54.6875 \text{ in}^2 \cdot 2 = 109.375$  [in<sup>2</sup>]

$37.13 \text{ in}^2 \cdot 2 = 75.46$  [in<sup>2</sup>]

Surface Area total =  $400.16 \text{ [in}^2] \cdot (1 \text{ [feet]}/12 \text{ [inches]})^2 = 25.8$  [feet<sup>2</sup>]


## Determine Irradiance Value

$141.6 \text{ [W]} / 25.8 \text{ [feet}^2] = 5.5 \text{ [W / feet}^2]$



Extreme case +40 degrees Fahrenheit  
Temperature with respect to ambient

- Conservative estimate
- Raspberry Pi has 178 degree Fahrenheit
- ANSYS Icepak for variable thermal loading
- No need to retrofit enclosure with cooling solutions
- Current analysis is undergoing final review from sponsoring professors



Cameron

# CONTROLLING THE PORTS

```
#NUTRIENT ANALYZER
#3/22/2022

import serial
import io
import time

ser = serial.Serial('/dev/ttyUSB0')
t0 = 0.001 # time for running pumps simultaneously
t1 = 2     # sleep time for valve movement
t2 = 8     # sytem flush dispense time (1000 ul at 150ul/sec)
t3 = 12    # time for pump movements with mixing (ex: 600 ul at 60ul/sec)
t4 = 16    # flow cell dispense time (400 ul at 25 ul/sec)
t5 = 300   # flow cell wait time for absorbance measurement (5 mins.)

# This section refers to port positions

def port_1():
    print("Moving to port 1 (Waste)")           # move to port 1 (Waste)
    ser.write(b'AMA 960\r\n')                  # Not used at this time
    time.sleep(t1)


def port_2():
    print("Moving to port 2 (Flow Cell)")       #Move to port 2 (flow cell)
    ser.write(b'AMA 720\r\n')
    time.sleep(t1)

def port_3():
    print("Moving to port 3 (Molybdate)")       # move to port 3 (Molybdate)
    ser.write(b'AMA 480\r\n')
    time.sleep(t1)

def port_4():
    print("Moving to port 4 (P04 Sample)")      # move to port 4 (P04 Sample)
    ser.write(b'AMA 240\r\n')
    time.sleep(t1)

def port_5():
    print("Moving to port 5 (Absorbic Acid Reagent)") #Move to port 5 (Absorbic Acid Reagent)
    ser.write(b'AMA 0\r\n')
    time.sleep(t1)

def port_6():
    print("Moving to port 6 (P04 Standard)")    # move to port 6 (P04 Standard)
    ser.write(b'AMA 1200\r\n')
    time.sleep(t1)
```





# CONTROLLING THE PUMPS

```
# This section is for system flush
def system_flush():
    print("System Flush In Progress")
    ser.write(b'CVM 150*EU\r\n')
    ser.write(b'CMR 1000*EU\r\n')
    time.sleep(t0)
    ser.write(b'DVM 150*EU\r\n')
    ser.write(b'DMR 1000*EU\r\n')
    time.sleep(t2)
    #This sets pump 1 to dispense 1000ul at 150 ul/sec
    #This sets pump 2 to dispense 1000ul at 150 ul/sec

#This section is for Dispenses
def blank_sample():
    print("Dispensing Blank Sample")
    ser.write(b'CVM 60*EU\r\n')
    ser.write(b'CMR 600*EU\r\n')
    time.sleep(t3)
    #This sets pump 1 to dispense 600ul at 60 ul/sec


def molybdate_reagent():
    print("Dispensing 320 ul pump_1/ aspirating 400 ul pump_2")
    ser.write(b'CVM 32*EU\r\n')
    ser.write(b'CMR 320*EU\r\n')
    time.sleep(t0)
    ser.write(b'DVM 40*EU\r\n')
    ser.write(b'DMR -400*EU\r\n')
    time.sleep(t3)
    #This sets pump 1 to dispense 320ul @ 32ul/sec
    #This sets pump 2 to aspirate 400 @ 40ul/sec

def ascorbic_acid_reagent():
    print("Aspirating 400 ul pump_1/ Dispensing 320 ul pump_2")
    ser.write(b'CVM 40*EU\r\n')
    ser.write(b'CMR -400*EU\r\n')
    time.sleep(t0)
    ser.write(b'DVM 32*EU\r\n')
    ser.write(b'DMR 320*EU\r\n')
    time.sleep(t3)
    #This sets pump 1 to aspirate 400 @ 40ul/sec
    #This sets pump 2 to dispense 320ul @ 32ul/sec

def flow_cell():
    print("Dispensing 400 ul pump_1 into flow cell")
    ser.write(b'CVM 25*EU\r\n')
    ser.write(b'CMR 400*EU\r\n')
    time.sleep(t4)
    time.sleep(t5)
    #This sets pump 1 to dispense 400ul at 25ul/sec

def po4_standard():
    print("Aspirating 600 ul pump_1 ")
    ser.write(b'CVM 60*EU\r\n')
    ser.write(b'CMR -600*EU\r\n')
    time.sleep(t3)
    #This sets pump 1 to aspirate 600ul at 60 ul/sec

def po4_sample():
    print("Aspirating 600 ul pump_1 ")
    ser.write(b'CVM 60*EU\r\n')
    ser.write(b'CMR -600*EU\r\n')
    time.sleep(t3)
    #This sets pump 1 to aspirate 600ul at 60 ul/sec
```



Tracy Hunter

# SEQUENCE OF CONTROL

```
#The system will run the blank sample sequence twice, then it will follow with the PO4 standard twice, then
#there will be one PO4 seawater sample run only once. The system will then post process to calculate the
#the PO4 concentration. The system will then sleep for ~40 mins., then run only the PO4 seawater sample again, followed by
#and additional sleep cycle. This will be repeated again, then calibration process will be restarted(Running the complete
#code from the beginning)

#(Maybe set up a for loop with sequences to follow specific sampling steps which will reset itself after 3 hours)
# ACTIVATE INSTRUMENT
#turn on analyzer
#turn off light
#get dark scan
#turn lamp on
#wait 3 minutes for lamp to warm up

#RUN BLANK SAMPLE (This will be done twice)
port_2()
system_flush()
#NEED REFERENCE SCAN HERE
blank_sample()
port_3()
molybdate_reagent()
port_5()
ascorbic_acid_reagent()
port_2()
flow_cell()
#ACQUIRE ABSORBANCE VALUE HERE

#RUN PO4 STANDARD (This will be done twice)
system_flush()
#NEED REFERENCE SCAN HERE
port_6()
po4_standard()
port_3()
molybdate_reagent()
port_5()
ascorbic_acid_reagent()
port_2()
flow_cell()
#ACQUIRE ABSORBANCE VALUE HERE

#RUN PO4 SAMPLE (This will be done once)
system_flush()
#NEED REFERENCE SCAN HERE
# RUN AUXILIARY PUMP HERE FOR 60 secs
port_4()
po4_sample()
port_3()
molybdate_reagent()
port_5()
ascorbic_acid_reagent()
port_2()
flow_cell()
#ACQUIRE ABSORBANCE VALUE HERE


#POST PROCESS HERE
ser.close()
```

This sequence is for one cycle only. Currently working to implement continuous application.

On 4/4/2022 there was an additional change to sequence by MLML.

Tracy Hunter

# Reagent Storage Design



## Reagent storage specs changed

**Initially:** Reagents would be hung from a bar separate from the electronics housing

**Scope change:** reagents need to be in the dark and submerged in water

## **Internal Dimension needed**

Width - 18"

Depth - 12"

Height - 22"

Tracy Hunter

# Container Solution

Home / Tamco® Industries / Tamco® Square & Rectangular Tanks / Hgt.




TAMCO INDUSTRIES

18" L x 12" W Black Standard  
Tamco® Tank Cover  
\$33.29



Shop Now



Home / Valves / Spigots & Fittings / Water Butt Tap

### Water Butt Tap

This water butt tap is durable and smooth, making it ideal for water cisterns, composting, horsebarns, and winemaking kits. It is also suitable for many home and garden applications.

- HDPE body & polypropylene spigot
- Gasket provides leak-proof seal (17351 only)
- Designed to fit through 1" hole in pail or drum sidewall
- High flow rate
- Handle turns independently of nozzle
- Smooth turning
- Meets FDA standards

Print


Stock	Description	Catalog Pg.	Price
17351	3/4" Black Long Thread Water Butt Tap with Gasket Save 5%: 25+   Save 10%: 50+   Save 15%: 150+	P-339	\$2.27/each In Stock
	3/4" Black Nut for Water Butt Tap Save 5%: 25+   Save 10%: 50+   Save 15%: 150+	P-339	\$0.96/each In Stock

[Water Butt Tap | U.S. Plastic Corp. \(usplastic.com\)](https://www.usplastic.com)

[20 Gallon Black Polyethylene Tank - 18" L x 12" W x 24" Hgt. | U.S. Plastic Corp. \(usplastic.com\)](https://www.usplastic.com)

Tracy Hunter







Item	Price
.375 X 48" X 96" HDPE SHEET	344.35
.375 X 54" X 96" Seaboard Marine Grade Sheet	404.32
375 X 48" X 96" ABS Sheet	241.68


Item	Price
20 Gallon Tank	120.67
Lid	33.29
Water Tap	2.67
Nut & Washer	0.97
Total	157.60


Tracy Hunter

# Enclosure Design



Andrew






Andrew










SOLIDWORKS Educational Product. For Instructional Use Only.


Andrew




# Enclosure Design




Andrew






Andrew





Andrew



SOLIDWORKS Educational Product. For Instructional Use Only.

Andrew



# Spectrometer Updates

- Initial issues
  - missing libraries
  - fixing file pathing
  - detecting hardware
- Software successfully installed and tested



Ocean Optic Flame S

Tim Ly



# Spectrometer Code

- Main functions Provided by Marine @MLML
- Our job:
  - split
  - integrate
  - test with lab values

```
#import necessary packages
import seabreeze
import pandas as pd
import numpy as np
import time
import matplotlib.pyplot as plt

#function for collecting dark scan
def darkscan(spec):
    """
    Inputs:
    -Spectrophotometer read from Seabreeze Library

    Outputs:
    -Intensity values across all wavelengths taken with lamp off.
    """

    time.sleep(10) #sleep 10s to ensure no light transmission after lamp turned off

    wavelengths = spec.wavelengths()
    intensities = spec.intensities()
    column_names = ['wavelengths', 'intensities']
    combine = np.vstack((wavelengths, intensities)).T
    np.shape(combine)


    darkscan = pd.DataFrame(data=combine, columns=column_names) #produces dataframe from most recent scan taken from spec.
```

Code provided by Marine from MLML



```
pi@0
File Edit Tabs Help
pi@0NA:~ $ python3 test.py
<Spectrometer FLAMES:FLMS02410>
```

Figure 1



## Next Steps for data collection

- Install the rest of the components
  - light source, flow cell, tubing
- Run tests using deionized water and blue dye
  - Compare numbers with values ran from lab deployed system
  - Adjust our system as needed





diagram of system flow

Parts subtasks		September					October					November					December					January					February					March					April					May				
		5	12	19	26		3	10	17	24	31	7	14	21	28		5	12	19	26		2	9	16	23	30	6	13	20	27		6	13	20	27		3	10	17	24		1	8	15	22	
Design / Planning																																														
Connecting	control of pumps																																													
	control of servo valve																																													
	control of spectrometer																																													
Hardware	User Interfacing, Terminal																																													
	Housing, Case																																													
System	combined controls																																													
	feedback method																																													
Implementation	Interfacing with existing facilities																																													
	Reagents / DI containment																																													
	Monitoring system for Supply & Power																																													
Interfacing	accomidations for wireless communications																																													
	display, controls, fuses, etc																																													
Manufacturing																																														
MCU																																														
Wiring																																														
process piping																																														
Programming																																														
Feedback loops																																														
Assembly																																														
soldering																																														
hardware																																														
process piping																																														
Reagent containment																																														
User Interface																																														
Validation																																														
Accuracy of Analysis																																														
continuous use test																																														
error flagging																																														
feedback tuning																																														
ancilaries																																														



SJSU

